



Issue date: 01/02/2026



Cod. 1045B

Bottle brush ø 50 mm – Blue

Bottle brush for cleaning of test tubes, bottles, tubes and tools with a cylindrical shape in general. Equipped with flexible handle for deep cleaning and ergonomic grip for a safe and comfortable hold. The anti-drip hole allows an optimal drying preserving the integrity of the bristle.



Certifications and regulatory compliance

- Regulation (EC) No 1907/2006 (REACH) and updates
- FDA CFR Title 21 and amendments
- Regulation (EU) No 1935/2004 and amendments
- Regulation (EU) No 10/2011 and amendments
- Regulation (EU) No 2023/2006 (Good Manufacturing Practices GMP) and amendments

Conditions of use and safety requirements¹

Operating temperature	min -20 / max 80 C°
Dishwasher resistance ²	Tested at 65°C according to BS EN 12875-1:2005
Steam sterilization	-
Instructions for use	• Remove all packaging and/or labels. • Sanitize/disinfect according to HACCP protocol. • Hang the tool using appropriate hooks/supports. • Allow to dry completely before use.
Signs of deterioration	The product must be regularly inspected. Observable signs of deterioration include: • Deformation of bristles (e.g. from dishwasher use); • Material softening, breakage, or discoloration; • Detachment of parts or structural damage. If any such issues are detected, the product must be withdrawn from food contact use and replaced, as also required under Regulation (EU) 2023/988 on general product safety.

Technical specifications

Length (mm)	500	Diameter (mm)	50
Width (mm)	-	Material	PP - SS - PBT
Height (mm)	-	Net weight (kg)	0,045
Bristle (mm)	Medium - Ø 0,5	Quantity per pack (pcs)	24
Bristle height (mm)	25		

¹ in accordance with Regulations (EU) 2025/351 and 2023/988

² Exposure to high temperatures or intensive cycles may cause bristle deformation or material degradation.